

Date 2/7/2013

Objective PCR amplification of triplet repeat region of FMR1 with a modified program (as below):
 - Without hot start
 - Three step amplification process
 - Denat temp grad 75-95deg
 - Extension and annealing temp.s and times (that have worked previously
 ie 54deg/ 1:30mins annealing-72deg/ 2:30mins extension)

Description *Template used:* 100ng genomic DNA NA20230, (53/54 rpt sample from Coriell)
 This is the newly ordered (Jan 2013) DNA: all old DNA is over FT for and FT rev
Primers: Herculaase Fusion II (Agilent)
Polymerase: Herculaase Fusion II (Agilent)
Solvents: No Solvt rxns
 Rxns in 0.5M NMP + 2.2M Betaine
 Rxns in 2.2M Betaine

<i>Reaction Composition:</i>		No solvt rxns	0.5M NMP + 2.2M Betaine rxns	2.2M Betaine rxns
Component	Final Conc/amt	Final Conc/amt	Final Conc/amt	
Herc Reac Buffer	1 X	1 X	1 X	
dNTPs	0.4 mM	0.4 mM	0.4 mM	
FT-For/ SFS-for	0.4 uM	0.4 uM	0.4 uM	
FT-Rev/ SFS-rev	0.4 uM	0.4 uM	0.4 uM	
gDNA NA203230	100 ng	100 ng	100 ng	
Taq Polymerase	0.25 ul	0.25 ul	0.25 ul	
Water				
NMP		0.5 M		
Betaine		2.2 M	2.2 M	
total	25ul	25ul	25ul	

Cycling Conditions:

- 1 75-95degC/ 30secs
- 2 54degC/ 1:30mins
- 3 72degC/ 2:30mins
- 4 GOTO 2 40 times


Gradient steps: 95, 94, 91.7, 87.7, 82.9, 79.1, 76.4, 75degC

Gel Electrophoresis: 10% precast polyacrylamide gels (Life Technologies)
 8ul of PCR rxn + 2ul of (5X) sample buffer loaded
 Molecular Ladders: 5ul loaded
 100bp ladder (NEB)
 50bp ladder (NEB)

Gel Pictures

Lanes are marked according to the Denat temp of the particular rxn
Expected product size is 383bp.

FT-No Solvt 20sec exposure



FT-Herc-0.5M NMP +2.2M Betaine 20sec exposure



FT-Herc + 2.2M Betaine 20sec exposure



- Comments** Expected bands not seen in any condition.
Even hot start rxns do not seem to have worked.
All Herc rxns were carried out with the new batch of gDNA.
(This expt was with the old (Ram's primer stks)