


Signal Test in Labelfree

31.8.16, ENC



signal:	1	2	3	4	5	6	7
sample:	buffer only	NAD+	Sirt3	Ac-peptide	DeAc-peptide	Honokiol	Honokiol
mM:		5	0.002	4.7	4.7	1	1
capillary type:	standard						premium coated zero background
buffer:	HDAC buffer, 0.1% Pluronic F-127			50% DMSO + HDAC buffer, 0.1% Pluronic F-127	HDAC buffer, 0.1% Pluronic F-127	10% DMSO + HDAC buffer, 0.1% Pluronic F-127	
%LED	20						